

Projeto MedVault

Equipe

- Líder

Jorge Abrão Neto

- Membros

David Samuel Tavares de Moraes

Lucas Mendes da Silva


Paulo Victor Rocha de Almeida



Objetivo

Objetivo inicial

Nosso projeto tinha como objetivo inicial, estudar e desenvolver um proxy de re-criptação para a distribuição de documentos médicos ligados ao armazenamento off chain acoplado a um sistema blockchain.




Objetivos

- Desenvolver uma aplicação para armazenamento seguro e descentralizado de registros médicos usando Blockchain e IPFS. ✓

- Desenvolvimento do proxy de re-criptação. ✓

- Desenvolvimento de contratos inteligentes para blockchain. ✓

- Integrar proxy à blockchain. ✗



Estrutura

IPFS



Solução escolhida para o armazenamento distribuído de dados. O IPFS é uma rede peer-to-peer que utiliza hashes criptográficos para endereçar conteúdo, permitindo que os dados sejam armazenados de forma descentralizada, recuperados eficientemente e resistam à censura. Essa arquitetura distribuída e resistente à censura não apenas melhora a velocidade e a disponibilidade do acesso aos dados, mas também reduz a carga sobre os pontos centrais, superando os desafios tradicionais de escalabilidade e centralização de armazenamento de dados.

```
49 def write(cid: str, path: str) -> None:
50     # Definindo os parâmetros para a solicitação GET no IPFS
51     params = {'arg': cid, 'archive': True, 'compress': True}
52
53     # Fazendo uma solicitação POST para obter um arquivo do IPFS
54     response = requests.post(f"{IPFS_BASE_URL}/get", params=params)
55
56     # Verificando se a resposta do IPFS foi bem-sucedida
57     if response.status_code != 200:
58         print(f"Erro ao obter arquivo do IPFS: {response.content}")
59         return
60
61     # Gravando o conteúdo da resposta em um arquivo local
62     with open(f"{path}/output.tar.gz", 'wb') as file:
63         file.write(response.content)
```

Elementos de robustez presentes



- **Tolerância**

O Hyperledger é projetado para ser tolerante a falhas, tanto de nós quanto de redes. Os nós do Sawtooth possuem mecanismos de detecção de falhas e são capazes de se recuperar e se reconfigurar em caso de falhas. Além disso, ele possui recursos para lidar com nós maliciosos ou desonestos, garantindo a segurança e a confiabilidade do sistema.

- **Coordenação**

A coordenação em um sistema Sawtooth é facilitada pelo mecanismo de validação de transações distribuídas. **Cada transação é submetida para validação em diferentes nós, e um algoritmo de consenso é usado para garantir que os nós cheguem a um acordo sobre a validade e a ordem das transações.** Isso requer coordenação eficiente entre nós para sincronizar e validar as transações de forma consistente.



Elementos de robustez presentes

- **Consenso**

O Sawtooth utiliza um mecanismo de consenso modular que significa que diferentes algoritmos de consenso podem ser implementados e escolhidos de acordo com as necessidades do sistema. O algoritmo de consenso Proof of Elapsed Time (PoET) foi implementado através do Sawtooth para, juntamente com alguns mecanismos de validação alcançar um consenso confiável e resistente a falhas.

- **Consistência**

A consistência dos dados é garantida no Sawtooth por meio de um modelo de transações em estado de finalização (Transaction Family State). Cada transação atualiza o estado do blockchain de maneira atômica, garantindo que as operações ocorram de forma consistente e correta. Além disso, o Proof Of Elapsed Time trabalha juntamente com com alguns mecanismos de validação, alcançar um consenso confiável e resistente a falhas.

Contratos Inteligentes

Smart Contracts:




- **O que são?**

Programas computacionais que executam automaticamente as regras acordadas previamente entre as partes. Sua lógica de funcionamento é a mesma de um documento legal tradicional. Ou seja, são estabelecidas obrigações, penalidades, benefícios e outras consequências devidas entre as pessoas que celebram o acordo.

- **Agilidade:** os contratos inteligentes são baseados em códigos de software, **Precisão:** sem preenchimentos manuais, as informações se tornam mais precisas. Além disso, elas são atualizadas automaticamente conforme o andamento das transações;
- **Autonomia:** por falar na ausência de terceiros e na automatização, isso elimina riscos e possíveis erros oriundos da manipulação por outros indivíduos;
- **Segurança:** os contratos inteligentes são criptografados e distribuídos pelos nós da blockchain. Isso elimina qualquer chance de perda ou alteração não autorizada;
- **Confiabilidade:** graças à autonomia e segurança dos smart contracts, suas operações são todas rastreáveis e há garantia de que as informações são legítimas.

Diagrama



Family_Controller

```
1 from sawtoothsdk.processor.handler import TransactionHandler
2
3 from models.doctor_model import Doctor
4 from models.patient_model import Patient
5 from utils import _hash
6
7 import json
8
9 FAMILY_NAME = 'FAMILY_CONTROLLER'
10
11 class ControllerTransactionHandler(TransactionHandler):
12     def __init__(self):
13         self._namespace_prefix = _hash(FAMILY_NAME.encode('utf-8'))[0:6]
14     @property
15     def family_name(self):
16         return FAMILY_NAME
17
18     @property
19     def family_versions(self):
20         return ['1.0']
21
22     @property
23     def namespaces(self):
24         return [self._namespace_prefix]
25
26     def getNamespace(self):
27         return self._namespace_prefix
28
29     def apply(self, transaction, context):
30         try:
31             header = transaction.header
32             signer = header.signer_public_key
33             action, payload = ControllerFactory.from_bytes(transaction.payload)
34             address = self._namespace_prefix + _hash(payload.cpf.encode('utf-8'))[:64]
35             state = context.get_state([address])
36         except:
37             print("Um erro ocorreu")
38             return
```



```
40     if action == 'add':
41         if state:
42             print(f'{payload._type} already exists')
43             return None
44         state_data = payload.to_bytes()
45         context.set_state({address: state_data})
46         print(f'{payload._type} was added')
47
48     elif action == 'show':
49         if not state:
50             print(f'{payload._type} does not exist')
51             return None
52         state_data = state[0].data.decode('utf-8')
53         print(f'{payload._type} data: {ControllerFactory.getPatient(state)}')
54
55     elif action == 'delete':
56         if not state:
57             print('{} does not exist' %format(payload._type))
58             return None
59         context.delete_state([address])
60         print(f'{payload._type} was deleted')
61
```



```
63 ~ class ControllerFactory:
64 ~     @staticmethod
65 ~     def getPatient(state):
66 ~         if not state:
67 ~             print('Patient does not exist')
68 ~             return None
69 ~         state_data = state[0].data
70 ~         return Patient.from_bytes(state_data)
71 ~     @staticmethod
72 ~     def getPayload(payload):
73 ~         try:
74 ~             # The payload is csv utf-8 encoded string
75 ~             data = json.loads(payload.decode())
76 ~             action = data["action"]
77 ~             type = data["type"]
78 ~             body = data["body"]
79 ~         except ValueError:
80 ~             print("Invalid payload serialization")
81 ~             return None
82 ~
83 ~         if not action:
84 ~             print('Action is required')
85 ~             return None
86 ~
87 ~         if action not in ('add', 'show', 'delete'):
88 ~             print('Invalid action: {}'.format(action))
89 ~             return None
90 ~
91 ~         if type not in ('doctor', 'patient'):
92 ~             print('Invalid type: {}'.format(type))
93 ~             return None
94 ~
95 ~         if type == 'doctor':
96 ~             return action, Doctor(body)
97 ~
98 ~         if type == 'patient':
99 ~             return action, Patient(body)
100 ~     @staticmethod
101 ~     def from_bytes(payload):
102 ~         return ControllerFactory.getPayload(payload=payload)
```

Record_Family

```

1  from sawtooth_sdk.processor.handler import TransactionHandler
2
3  from models.record_model import Record
4  from utils import _hash
5  from families.controller_family import ControllerTransactionHandler, ControllerFactory
6
7  import json
8
9  FAMILY_NAME = 'FAMILY_RECORD'
10
11 class RecordTransactionHandler(TransactionHandler):
12     def __init__(self):
13         self._namespace_prefix = _hash(FAMILY_NAME.encode('utf-8'))[0:6]
14
15     @property
16     def family_name(self):
17         return FAMILY_NAME
18
19     @property
20     def family_versions(self):
21         return ['1.0']
22
23     @property
24     def namespaces(self):
25         return [self._namespace_prefix]
26
27     def apply(self, transaction, context):
28         header = transaction.header
29         signer = header.signer_public_key
30         action, patient_cpf, body = RecordFactory.from_bytes(transaction.payload)
31
32         controller_namespace_prefix = ControllerTransactionHandler().getNamespace()
33         patient_address = controller_namespace_prefix + _hash(patient_cpf.encode('utf-8'))[:64]
34         state = context.get_state([patient_address])
35         patient = ControllerFactory.getPatient(state)
36

```





```
37     if action == 'add':
38         record = Record(body)
39         patient.add_record(record)
40         print(f"Record added")
41     elif action == 'show':
42         id_record = body["id_record"]
43         record = patient.get_record(id_record)
44         print(f"Record {repr(record)}")
45     elif action == 'delete':
46         id_record = body["id_record"]
47         patient.delete_record(id_record)
48     elif action == 'grant':
49         id_record = body["id_record"]
50         doctor_cpf = body["doctor_cpf"]
51         id_request = body["id_request"]
52         request_status = body["request_status"]
53         patient.grant_record(id_record, doctor_cpf, id_request, request_status)
54         print(f"Record granted")
55         record = patient.get_record(id_record)
56         print(f"Record {repr(record)}")
57     elif action == 'request':
58         id_record = body["id_record"]
59         doctor_cpf = body["doctor_cpf"]
60         id_request = body["id_request"]
61         patient.request_record(id_record, doctor_cpf, id_request)
62         print(f"Record requested")
63     patient.update_patient(state=state, address=patient_address, context=context)
64     print(patient)
```

```
66 class RecordFactory:
67     @staticmethod
68     def getPayload(payload):
69         try:
70             # The payload is csv utf-8 encoded string
71             data = json.loads(payload.decode())
72             action = data["action"]
73             patient_cpf = data["patient_cpf"]
74             body = data["body"]
75
76         except ValueError:
77             print("Invalid payload serialization")
78             return None
79
80         if not action:
81             print('Action is required')
82             return None
83
84         if action not in ('add', 'show', 'delete', 'grant', 'request'):
85             print(f'Invalid action: {action}')
86             return None
87         return action, patient_cpf, body
88     @staticmethod
89     def from_bytes(payload):
90         return RecordFactory.getPayload(payload=payload)
```



Entidades

```
1 import json
2 from utils import _hash
3
4 class Doctor:
5     def __init__(self, body):
6         cpf = body.get("cpf", None)
7         name = body.get("name", None)
8
9         if not cpf:
10             print('CPF is required')
11             return None
12
13         if not name:
14             print('Name is required')
15             return None
16
17         self._name = name
18         self._cpf = cpf
19         self._type = "Doctor"
20
21     @property
22     def cpf(self):
23         return self._cpf
24
25     @property
26     def name(self):
27         return self._name
28
29     def to_bytes(self):
30         doctor = {
31             "name": self._name,
32             "type": self._type,
33             "cpf": self._cpf
34         }
35
36         return json.dumps(doctor).encode()
```





```
1 import json
2 from utils import _hash
3
4 from models.record_model import Record
5
6 class Patient:
7     def __init__(self, body):
8         cpf = body.get("cpf", None)
9         name = body.get("name", None)
10        records = body.get("records", None)
11
12        if not cpf:
13            print('CPF is required')
14            return None
15
16        if not name:
17            print('Name is required')
18            return None
19
20        if records:
21            records = [Record.from_json(record) for record in records]
22
23        self._name = name
24        self._cpf = cpf
25        self._records = records or []
26        self._type = "Patient"
27
28        @property
29        def cpf(self):
30            return self._cpf
31
32        @property
33        def name(self):
34            return self._name
35
36        @staticmethod
37        def from_bytes(data):
38            body = json.loads(data.decode('utf-8'))
39            return Patient(body)
40
41        def to_bytes(self):
42            patient = {
43                "name": self._name,
44                "cpf": self._cpf,
45                "type": self._type,
46                "records": [record.to_json() for record in self._records]
47            }
48            return json.dumps(patient).encode()
```

```
50 def update_patient(self, state, address, context):
51     state_data = self.to_bytes()
52     context.set_state({address: state_data})
53
54 def add_record(self, record):
55     self._records.append(record)
56
57 def get_record(self, id_record):
58     for record in self._records:
59         if record._id_record == id_record:
60             return record
61     print('Record does not exist')
62     return None
63
64 def delete_record(self, id_record):
65     for record in self._records:
66         if record._id_record == id_record:
67             self._records.remove(record)
68             return None
69     print('Record does not exist')
70     return None
71
72
73 def request_record(self, id_record, doctor_cpf, id_request):
74     for record in self._records:
75         if record._id_record == id_record:
76             record.requested(id_request, doctor_cpf)
77             return None
78     print('Record does not exist')
79     return None
80
81 def grant_record(self, id_record, doctor_cpf, id_request, request_status):
82     for record in self._records:
83         if record._id_record == id_record:
84             record.granted(id_request, request_status)
85             return None
86     print('Record does not exist')
87     return None
88
89 def __repr__(self) -> str:
90     return f"Name: {self._name}, CPF: {self._cpf}, Type: {self._type}, Records: {';'.join([repr(record) for record in self._records])}\n"
```




```
1 import json
2 from utils import _hash
3 from models.request_model import Request
4 class Record:
5     def __init__(self, body):
6         id_record = body.get("id_record", None)
7         title = body.get("title", None)
8         bundle_hash = body.get("bundle_hash", None)
9         requests = body.get("requests", None)
10
11         if not id_record:
12             print('ID is required')
13             return None
14
15         if requests:
16             requests = [Request.from_json(request) for request in requests]
17
18         self._id_record = id_record
19         self._title = title
20         self._bundle_hash = bundle_hash
21         self._requests = requests or []
22     @property
23     def action(self):
24         return self._action
25
26     @property
27     def title(self):
28         return self._title
29
30     @property
31     def bundle_hash(self):
32         return self._bundle_hash
33
34     def to_json(self):
35         record = {
36             "id_record": self._id_record,
37             "title": self._title,
38             "bundle_hash": self._bundle_hash,
39             "requests": [request.to_json() for request in self._requests]
40         }
41
42         return json.dumps(record)
```



```
44 @staticmethod
45 def from_json(record):
46     return Record(json.loads(record))
47
48 def requested(self, id_request, doctor_cpf):
49
50     request = Request({
51         "request_id": id_request, "doctor_cpf": doctor_cpf
52     })
53     self._requests.append(request)
54
55 def granted(self, id_request, request_status):
56     for request in self._requests:
57         if request._request_id == id_request:
58             request.reply(request_status)
59             return None
60     print('Record does not exist')
61     return None
62 def __repr__(self):
63     return f"ID: {self._id_record}, bundle_hash: {self._bundle_hash}, title: {self._title}, Requests: --- {' '.join([repr(request) for request in self._requests])} ---\n"
```

```
2 ~ import json
3   from utils import _hash
4
5 ~ class Request:
6 ~     def __init__(self, body):
7         doctor_cpf = body.get("doctor_cpf", None)
8         request_id = body.get("request_id", None)
9         request_status = body.get("request_status", None)
10
11         if request_status == None:
12             request_status = 0
13
14         self._doctor_cpf = doctor_cpf
15         self._request_id = request_id
16         self._request_status = request_status #0 for waiting approval, 1 for approved, 2 for denied.
17
18     @property
19     def doctor(self):
20         return self._doctor_cpf
21
22     @property
23     def id(self):
24         return self._request_id
25
26     @property
27     def status(self):
28         return self._request_status
29
30     def reply(self, status):
31         self._request_status = status
32
33     def to_json(self):
34         request = {
35             "doctor_cpf": self._doctor_cpf,
36             "request_id": self._request_id,
37             "request_status": self._request_status
38         }
39
40         return json.dumps(request)
```





```
42 @staticmethod
43 def from_json(data):
44     data_dict = json.loads(data)
45
46     # Extract the values for doctor_cpf and request_id
47
48     request = Request(data_dict)
49
50     return request
51
52 def __repr__(self):
53     return f"ID: {self._request_id}, Doctor_CPF: {self._doctor_cpf}, request_status: {self._request_status}"
```

Re-criptografia

Interface IPFS método add:



```
5 def add(file) -> str:
6     # Faz uma solicitação POST para adicionar um arquivo ao IPFS
7     response = requests.post(f"{IPFS_BASE_URL}/add", files={"file": file}, params={'cid-version': 1})
8
9     # Verifica se a resposta do IPFS foi bem-sucedida
10    if response.status_code != 200:
11        print(f"Erro ao obter arquivo do IPFS: {response.content}")
12        return
13
14    # Analisa os dados JSON retornados pela resposta
15    data = json.loads(response.content)
16
17    # Retorna o hash do CID (Content Identifier) do arquivo adicionado ao IPFS
18    return data["Hash"]
```

Interface IPFS método get:



```
20 def get(cid: str) -> dict:
21     # Parâmetros para a solicitação GET no IPFS
22     params = {'arg': cid, 'archive': True, 'compress': True}
23
24     # Faz uma solicitação POST para obter um arquivo do IPFS
25     response = requests.post(f"{IPFS_BASE_URL}/get", params=params)
26
27     # Verifica se a resposta do IPFS foi bem-sucedida
28     if response.status_code != 200:
29         print(f"Erro ao obter arquivo do IPFS: {response.content}")
30         return
31
32     try:
33         # Importa as bibliotecas necessárias para manipulação de arquivos tar e io
34         import tarfile, io
35
36         # Abre o arquivo tar recebido como resposta e lê o seu conteúdo
37         with tarfile.open(fileobj=io.BytesIO(response.content), mode='r:gz') as tar:
38             filename = tar.getnames()[0]
39             file = tar.extractfile(filename)
40             content: bytes = file.read()
41
42             # Converte o conteúdo lido (bytes) para um dicionário usando JSON
43             return json.loads(content)
44     except tarfile.BadGzipFile:
45         print("Erro ao obter arquivo do IPFS: Não foi possível realizar leitura do arquivo - Arquivo não é Gzip")
46     except KeyError:
47         print(f"Erro ao obter arquivo do IPFS: Não foi possível extrair arquivo {filename} - Arquivo não existe")
```

Interface IPFS método write:



```
49 def write(cid: str, path: str) -> None:
50     # Definindo os parâmetros para a solicitação GET no IPFS
51     params = {'arg': cid, 'archive': True, 'compress': True}
52
53     # Fazendo uma solicitação POST para obter um arquivo do IPFS
54     response = requests.post(f"{IPFS_BASE_URL}/get", params=params)
55
56     # Verificando se a resposta do IPFS foi bem-sucedida
57     if response.status_code != 200:
58         print(f"Erro ao obter arquivo do IPFS: {response.content}")
59         return
60
61     # Gravando o conteúdo da resposta em um arquivo local
62     with open(f"{path}/output.tar.gz", 'wb') as file:
63         file.write(response.content)
```


Encriptação do token de acesso (simula o trabalho da blockchain):



```
16 def create_instance(content: bytes) -> None:
17     patient_keys = keygen()
18     doctor_keys = keygen()
19     patient_signer_keys = keygen()
20     patient_signer = Signer(patient_signer_keys['secret_key'])
21     doctor_id = "1"
22     record = "RECORD 1"
23     key_frags = patient_grants_access(patient_keys['secret_key'], patient_signer, doctor_keys['public_key'], 1, 1)
24
25     cid_record = patient_create_record(content, patient_keys['public_key'])
26     print(f"O Record do paciente foi armazenado na IPFS em {cid_record}")
27
28     data = {
29         'doctor_id': doctor_id,
30         'public_key_patient': serializeFrom(bytes(patient_keys['public_key'])),
31         'public_key_signer_patient': serializeFrom(bytes(patient_signer_keys['public_key'])),
32         'public_key_doctor': serializeFrom(bytes(doctor_keys['public_key'])),
33         'key_frag': serializelistFrom(key_frags),
34         'link': cid_record
35     }
36     # Armazenando o pedido do SC em um arquivo
37     with open(SC_FILE_STORAGE, 'w') as file:
38         json.dump(data, file, indent=4)
39     print("Arquivo do Contrato Inteligente foi criado")
40
41     # Criando o armazenamento do Doutor:
42     with open(DOCTOR_FILE_STORAGE, 'w') as file:
43         keys = serializeKeys(doctor_keys)
44         json.dump(keys, file, indent=4)
45     print("Arquivo com as chaves do Doutor foi criado")
46
```

Re-criptação:



```
47 def test_instance():
48     # Armazenando o pedido do SC em um arquivo
49     with open(SC_FILE_STORAGE, 'r') as file:
50         token_obj = json.load(file)
51     print("Token obtido do armazenamento do SC")
52
53     # Criando o armazenamento do Doutor:
54     with open(DOCTOR_FILE_STORAGE, 'r') as file:
55         doctor_key_encoded = json.load(file)
56         doctor_keys = deserializeKeys(doctor_key_encoded)
57     print("Chaves obtidas do armazenamento do doutor")
58
59     public_key_patient = PublicKey.from_bytes(deserializeFrom(token_obj["public_key_patient"]))
60     public_key_signer_patient = PublicKey.from_bytes(deserializeFrom(token_obj["public_key_signer_patient"]))
61     public_key_doctor = PublicKey.from_bytes(deserializeFrom(token_obj["public_key_doctor"]))
62     key_frag = deserializeListFrom(VerifiedKeyFrag.from_verified_bytes, token_obj["key_frag"])
63     print("Token decodificado")
64
65
66     ipfs_obj = ipfs.get(token_obj["link"])
67     capsule = Capsule.from_bytes(deserializeFrom(ipfs_obj['capsule']))
68     ciphertext = deserializeFrom(ipfs_obj['content'])
69     print("Cifra e capsula obtidas do IPFS")
70
71     cfrags = get_cfrags(public_key_signer_patient, public_key_patient, public_key_doctor, capsule, kfrags=key_frag)
72     content: bytes = decrypt_reencrypted(
73         verified_cfrags = cfrags,
74         receiving_sk = doctor_keys["secret_key"], delegating_pk = public_key_patient,
75         capsule = capsule, ciphertext = ciphertext
76     )
77     print("Conteudo decodificado por recriptacao")
78     return content
```

Demonstração



Docker Desktop

Update to latest

Search for local and remote images, containers, and more...

Ctrl+K

Sign in

Containers

Images

Volumes

Dev Environments BETA

Docker Scout EARLY ACCESS

Learning Center

Extensions

Add Extensions

Containers [Give feedback](#)

Your running containers show up here

A container is an isolated environment for your code, run a container to see it here

What is a container?

Hands-on guide (5 mins)

Run Docker Hub images

Hands-on guide (5 mins)

```
1 FROM node
2 RUN mkdir -p
3 WORKDIR /app
4 COPY package
```

How do I run a container?

Hands-on guide (6 mins)

Publish your image

Hands-on guide (5 mins)

[More samples and guides at the Learning Center](#)

Not connected to Hub

v4.20.1



```
PS C:\sawtooth-example\sawtooth> docker compose -f sawtooth-default.yaml up -d
```

```
[+] Building 0.0s (0/0)
```

```
[+] Running 8/8
```

✓ Container sawtooth-validator-default	Started	0.5s
✓ Container sawtooth-devmode-engine-rust-default	Started	2.1s
✓ Container sawtooth-intkey-tp-python-default	Started	2.0s
✓ Container sawtooth-settings-tp-default	Started	2.1s
✓ Container sawtooth-rest-api-default	Started	1.5s
✓ Container sawtooth-xo-tp-python-default	Started	1.5s
✓ Container sawtooth-shell-default	Started	2.6s
✓ Container sawtooth-dev-default	Started	2.4s

```
PS C:\sawtooth-example\sawtooth> █
```

```
PS C:\sawtooth-example\sawtooth> docker exec -t sawtooth-dev-default bash -c "cd sawtooth/server/ ; python3 main.py"
```

```
Servidor iniciou
```

```
█
```



Containers



Images



Volumes

Dev Environments BETADocker Scout EARLY ACCESS

Learning Center

Extensions



Add Extensions

Containers

[Give feedback](#)

Container CPU usage ⓘ

0.66% / 800% (8 cores allocated)

Container memory usage ⓘ

122.06MB / 5.6GB

Show charts ▾



Search



Only show running containers



Name

Image

Sta... ↑



Port(s)

Last started

CPU (%)

Actions

[sawtooth](#)

Running (6/8)

5 minutes ago

0.66%



Showing 1 item

```

rawPayload = {
    "action": "add",
    "type": "doctor",
    "body": {
        "cpf": DOCTOR_CPF,
        "name": "ryan",
    }
}

client = PIBITIClient(baseUrl=DEFAULT_URL, keyFile=key_file, family_name=CONTROLLER_FAMILY_NAME, cpf=DOCTOR_CPF)
print("Wrap and send - Add Doctor")
response = client._wrap_and_send(rawPayload)
print("Response: {}".format(response))

```

```
PS C:\sawtooth-example\sawtooth> docker exec -t sawtooth-dev-default bash -c "cd sawtooth/client/ ; python3 dapp_test.py"
```

```
Wrap and send - Add Doctor
```

```

Response: {
  "link": "http://rest-api:8008/batch_statuses?id=2befb43d15d6cf44e8004c12be2ece209ed486b8a2ba081f2752a888a7defb08549e75fb49f285a92db79bed251d664156d821d14411666cb0860867106af821"
}

```

```
PS C:\sawtooth-example\sawtooth> █
```

```
docker exec -t sawtooth-dev-default bash -c "cd sawtooth/server/ ; python3 main.py"
```

```
Servidor iniciou
```

```
Doctor was added
```

```
Doctor was added
```

```
█
```

```

rawPayload = {
    "action": "add",
    "type": "patient",
    "body": {
        "cpf": PATIENT_CPF,
        "name": "Joana",
    }
}

client = PIBITIClient(baseUrl=DEFAULT_URL, keyFile=key_file, family_name=CONTROLLER_FAMILY_NAME, cpf=PATIENT_CPF)
print("Wrap and send - Add Patient")
response = client._wrap_and_send(rawPayload)
print("Response: {}".format(response))

```

```

PS C:\sawtooth-example\sawtooth> docker exec -t sawtooth-dev-default bash -c "cd sawtooth/client/ ; python3 dapp_test.py"

```

```

Wrap and send - Add Patient

```

```

Response: {
  "link": "http://rest-api:8008/batch_statuses?id=b1047d02d84bc622a571af937ef81d0469239d077c27177856faedd70bfa4fdf4df6fc0f1ed5b94792c129f5bd2f7c89b1272d8d0823a0454e4282d10d250f24"
}
PS C:\sawtooth-example\sawtooth> 

```

```

docker exec -t sawtooth-dev-default bash -c "cd sawtooth/server/ ; python3 main.py"

```

```

Servidor iniciou
Doctor was added
Doctor was added
Patient was added
Patient was added

```



```

rawPayload = {
    "action": "add",
    "patient_cpf": PATIENT_CPF,
    "body": {
        "id_record": "1",
        "title": "Record Title 1",
        "bundle_hash": "BUNDLE_HASH12312312",
    }
}

client = PIBITIClient(baseUrl=DEFAULT_URL, keyFile=key_file, family_name=RECORD_FAMILY_NAME, cpf=PATIENT_CPF)
print("Wrap and send - Add Record")
response = client._wrap_and_send(rawPayload)
print("Response: {}".format(response))

```

```

PS C:\sawtooth-example\sawtooth> docker exec -t sawtooth-dev-default bash -c "cd sawtooth/client/ ; python3 dapp_test.py"
Wrap and send - Add Record
Response: {
  "link": "http://rest-api:8008/batch_statuses?id=886a9642ba003965caa9944752c0426bec72d8786eda75b1286de77d2d34f4c60f56ba6b63a1d3f7d61910f41a68db45ac93ae629a25ea33c9dcdd90587ebdd"
}
PS C:\sawtooth-example\sawtooth>

```

```

Doctor was added
Doctor was added
Patient was added
Patient was added
Record added
Name: Joana, CPF: 222.222.222-22, Type: Patient, Records: ID: 1, bundle_hash: BUNDLE_HASH12312312, title: Record Title 1, Requests: --- ---

Record added
Name: Joana, CPF: 222.222.222-22, Type: Patient, Records: ID: 1, bundle_hash: BUNDLE_HASH12312312, title: Record Title 1, Requests: --- ---

```

```

rawPayload = {
    "action": "request",
    "patient_cpf": PATIENT_CPF,
    "body": {
        "id_record": "1",
        "doctor_cpf": DOCTOR_CPF,
        "id_request": "ID_rREQUEST1",
    }
}

client = PIBITIClient(baseUrl=DEFAULT_URL, keyFile=key_file, family_name=RECORD_FAMILY_NAME, cpf=PATIENT_CPF)
print("Wrap and send - Request Record")
response = client._wrap_and_send(rawPayload)
print("Response: {}".format(response))

```

```

PS C:\sawtooth-example\sawtooth> docker exec -t sawtooth-dev-default bash -c "cd sawtooth/client/ ; python3 dapp_test.py"

```

```

Wrap and send - Request Record

```

```

Response: {
  "link": "http://rest-api:8008/batch_statuses?id=a0d582614a96ac94cf787789380c4e1060e4b02eff05b0578090b5ece69888b82a7238601f682ded8d0681b69ca303344cf17f
ee1dbe34bd155113a67a140aa1"
}

```

```

PS C:\sawtooth-example\sawtooth> 

```

Record requested

Name: Joana, CPF: 222.222.222-22, Type: Patient, Records: ID: 1, bundle_hash: BUNDLE_HASH12312312, title: Record Title 1, Requests: --- ID: ID_rREQUEST1, Doctor_CPF: 111.111.111-11, request_status: 0 ---

Record requested

Name: Joana, CPF: 222.222.222-22, Type: Patient, Records: ID: 1, bundle_hash: BUNDLE_HASH12312312, title: Record Title 1, Requests: --- ID: ID_rREQUEST1, Doctor_CPF: 111.111.111-11, request_status: 0 ---



IPFS



STATUS



FILES



EXPLORE



PEERS



SETTINGS

UI v2.15.0
Revision fadd900
[See the code](#)
[Report a bug](#)

QmHash/bafyHash

Browse



Files

4B
FILES

7MiB
ALL BLOCKS

+ Import

Name ↑


Pin Status

Size


No files in this directory. Click the "Import" button to add some.

No files here yet! Add files to your local IPFS node by clicking the **Import** button above.





EPSON Easy Photo Print Photo Print





```
81 def main():
82     with open("./examples/Hello World.pdf", "rb") as file:
83         content = file.read()
84         create_instance(content)
85
86     content = test_instance()
87     with open("./examples/pdf-example.pdf", "wb") as file:
88         file.write(content)
89     print("Conteudo escrito em um arquivo")
90
91 if(__name__=="__main__"):
92     main()
```

```
PS C:\Users\jorge\OneDrive\Área de Trabalho\codes\reencryption-proxy-master> python setup.py
O Record do paciente foi armazenado na IPFS em QmYfQzioQL5RP7MP5tWZz48YV2on8545YaDHcy7G6XNHAc
```

Arquivo do Contrato Inteligente foi criado

Arquivo com as chaves do Doutor foi criado

Token obtido do armazenamento do SC



Token decodificado


Cifra e capsula obtidas do IPFS

Conteudo decodificado por recriptacao

Conteudo escrito em um arquivo

```
PS C:\Users\jorge\OneDrive\Área de Trabalho\codes\reencryption-proxy-master> █
```







IPFS



STATUS



FILES



EXPLORE



PEERS



SETTINGS

UI v2.15.0

Revision fadd900

[See the code](#)

[Report a bug](#)

QmHash/bafyHash

Browse



Files

15KiB
FILES

8MiB
ALL BLOCKS

+ Import

Name ↑

Pin Status

Size

 Hello World.pdf
QmYFQzi0QL5RP7MP5tMz48YV2on8545YaDHcy7G6XNIHAc



15 KiB





IPFS



FILES



EXPLORE



SETTINGS



Files / Hello World.pdf

15KiB
FILES

8MiB
ALL BLOCKS

... More



Hello World

1 / 1

/ 1

100%



1

Hello World!
Exemplo de exame
aqui!-----

Referências

Referências Bibliográficas

Aplicação para armazenamento seguro e descentralizado de registros médicos usando Blockchain e IPFS

- [1] Madine, M. M., Al-Ayyoub, M., Jararweh, Y., & Al-Zoubi, H. (2020). Blockchain for Giving Patients Control Over Their Medical Records. IEEE Access, 8, 217487-217500.
- [2] Kumar, S., Bharti, A. K., & Amin, R. (2021). Decentralized secure storage of medical records using Blockchain and IPFS: A comparative analysis with future directions. Security and Privacy, 4(5), e173.
- [3] Kumar, S., Bharti, A. K., & Amin, R. (2021). Secure and Scalable Decentralized Supply Chain Management Using Blockchain and IPFS. IEEE Communications Magazine, 59(10), 68-74.



Obrigado!